

.

---

**From:** Piro, Peter (DPH)  
**Sent:** Monday, June 07, 2010 1:10 PM  
**To:** Hanchett, James (DPH)  
**Subject:** RE: TFMPP

Hi Jim,  
Do you use a 0.25 ID column or 0.32 ID? We seem to get peak splitting for the salt form unless we convert it to the base.  
Thanks  
Peter

---

**From:** Hanchett, James (DPH)  
**Sent:** Monday, April 26, 2010 8:16 AM  
**To:** Piro, Peter (DPH)  
**Subject:** TFMPP

Hi Peter,  
Talked with Sonja, she had TFMPP (Triflouromehtylphenyl piperazine) in some of her MDMA's. HP  
1MS gave the best separation about 12 seconds. GC condition are 120 degrees for 2.50 minutes then ramp @ 10  
degree/min to 280. MDMA elutes @ around 8 minutes, with TFMPP eluting right before. Try this site for info.  
<http://en.wikipedia.org/wiki/Trifluoromethylphenylpiperazine>

Jim